http://www.lenstec.com  
ISO13485 Registered Company


**LENSTEC (BARBADOS) INC.**, Airport Commercial Centre, Pilgrim Road, Christ Church BB17092, Barbados  
 ☎ +1 (246) 420-6795 📠 +1 (246) 420-6797 ✉ lenstecbarbados@lenstec.com

**EC REP** CMI spol. s.r.o. Trenčianska 47,821 09 Bratislava, Slovakia  
 ☎ +421 (2) 52631441 📠 +421 (2) 52922882 ✉ lenstec@cmi.sk

**UK SALES: LENSTEC (BARBADOS) INC.**, Lenstec House, Unit 8, Mariner Court, Calder Park, Wakefield, WF4 3FL, UK  
 ☎ +44 (0) 1924 382 678 📠 +44 (0) 1924 850 454 ✉ lenstecuk@lenstec.com


**WORLD HEADQUARTERS: LENSTEC, INC.**, 1765 Commerce Ave N, Saint Petersburg, Florida 33716, U.S.A.  
 ☎ +1 (727) 571-2272 📠 +1 (727) 571-1792 ✉ lenstec@lenstec.com

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**IMPORTANT NOTICE**

It is highly recommended that the surgeon adheres to the recommendations, precautions, contraindications and warnings outlined in these instructions.

**DEVICE DESCRIPTION**

The LENSTEC SOFTEC series of posterior chamber intraocular lenses (PCIOL) are ultraviolet absorbing, single-piece "C" loop intraocular lenses intended for the replacement of the human crystalline lens following phacoemulsification cataract removal. The LENSTEC SOFTEC intraocular lenses are manufactured from a medical grade co-polymer of hydrophilic acrylic, with a polymerisable UV blocker. The hydrophilic nature of the SOFTEC material reduces problems associated with silicone oil adhesion and silicone oil induced opacification<sup>2-4</sup>. The SOFTEC series of IOLs are designed with square edges<sup>5</sup>.

**INTENDED USE**

The LENSTEC SOFTEC series of posterior chamber lenses are intended for the replacement of the human crystalline lens following phacoemulsification cataract removal in adults 18 years and older. The lens is indicated for capsular bag placement.

**CONTRAINDICATIONS**

Outside of general contraindications for ocular surgery, the following specific contraindications apply: Uncontrolled glaucoma, microphthalmia, chronic severe uveitis, retinal detachment, corneal decompensation, diabetic retinopathy, iris atrophy, perioperative complications, potential foreseeable postoperative complications and other conditions which an ophthalmic surgeon might identify based on their experience.

**CLINICAL BENEFITS**

The primary benefits for clinical management and patient health include the treatment of:

- Aphakia
- Cataract
- Myopia
- Hyperopia

The duration of the treatment effect is anticipated to be permanent.

The Softec HD PLI produces the following clinical benefits:

- Reduces handling and loading errors
- Minimises the chance of infection
- Reduces operating time turnover

- Consistent, predictable and controlled insertion
- Prevention of loss of sight, and improved visual acuity

**PERFORMANCE CHARACTERISTICS**

The anticipated effects on patient quality of life are prevention of loss of sight and improved visual acuity.

**RISKS**

The potential risks of using the Lenstec Softec HD PLI are as follows:

Endophthalmitis, Toxic anterior segment syndrome (TASS), Anterior capsule fibrosis, Uveitis glaucoma hyphaema syndrome, Iritis, Iris capture, Cystoid macular edema, Corneal stromal edema, Posterior capsular contraction & lens deformation, Capsular damage, Decentration/tilt (small optic), Elevated IOP, Concomitant surgery, Implant material clouding

**WARNINGS**

The implanting ophthalmic surgeon shall consider the following warnings, and identify a risk/benefit ratio prior to surgery:

1. Failure to follow the implantation instructions supplied with this lens could lead to mishandling and subsequent IOL damage prior to or during implantation.
2. There is no clinical data to support placing this lens in the ciliary sulcus.
3. Any posterior capsulotomy opening should be limited to approximately 4 mm. Consistent with other IOLs, there is an increased risk of lens dislocation and/or secondary surgical reintervention with early or large YAG capsulotomies.
4. The SOFTEC series of intraocular lenses should not be implanted if the capsular bag is not intact or if there is significant zonular rupture/dehiscence.
5. The effectiveness of ultraviolet light absorbing lenses in reducing the incidence of retinal disorders has not been established. As a precaution, patients should be informed that they should wear sunglasses with UV protection when in sunlight.
6. The rate of cystoid macular edema may increase with extracapsular bag placement of the haptics.
7. Patients with any of the following could be at increased risk for complication(s) following implantation of any of the SOFTEC series of IOLs: previous ocular surgery, those meeting any of the listed factors in the 'Contraindications' section of this document, non-age related cataract, vitreous loss, iris atrophy, severe aniseikonia, ocular hemorrhage, macular degeneration or suspected microbial infection.
8. Patients who present complications at the time of cataract extraction could be at increased risk for complication(s) following implantation of any of the SOFTEC series of IOLs. This may include, but is not limited to: persistent bleeding, significant iris damage, uncontrolled positive pressure or significant vitreous prolapse or loss.
9. Patients Whether intraocular lens implantation would deleteriously affect the surgeon's ability to otherwise observe, diagnose or treat posterior segment diseases in the patient.
10. Whether patients who have a distorted eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible should have any of the SOFTEC series of IOLs implanted.
11. Whether patients who have recurrent severe anterior or posterior segment inflammation or uveitis should have any of the SOFTEC series of IOLs implanted.
12. Any circumstances which could lead to damage to the corneal endothelium during implantation should be avoided.
13. Children under the age of 2 are not suitable candidates for intraocular lenses.
14. Re-use of the IOL is strictly prohibited, as it raises serious safety and efficacy concerns:
  - LENSTEC does not provide cleaning/sterilization instructions. An improperly cleaned and/or sterilized IOL can cause significant damage to a patient's vision, due in part to cross

contamination induced infection.

- Once removed from its original packaging, the IOL can lose traceability. In the event an IOL is re-used, it is unlikely the user will know the correct expiry date, serial number or dioptric power.
- LENSTEC can not guarantee stability or proper function of either haptic or optic portions in the event that an IOL is re-used. Failure of either of these components can render the IOL ineffective.
- The injector is designed for single-use and components of the device are not able to be re-used. Attempted re-use of the injector will result in damage to the IOL which could cause serious harm to the patient.

#### **PRECAUTIONS**

1. Do not attempt to re-use the lens. Do not autoclave or attempt to re-sterilize the lens. Lenses requiring re-sterilization should be returned to LENSTEC.
2. Do not use the device if sterile packaging has been damaged or if there are traces of leakage on the pack.
3. Do not soak the intraocular lens with any solution other than a sterile balanced salt solution.
4. Once packaging has been opened, the intraocular lens must be used immediately. The hydrophilic nature of the lens can cause it to absorb substances with which it comes into contact, such as disinfectants, medicines, blood cells, etc. This may cause a "Toxic Lens Syndrome".
5. The lens must be implanted within 2 minutes following removal from its saline bath, as dehydration causes the lens material to become brittle.
6. The lens must be implanted in the capsular bag.
7. Do not use the intraocular lens after the expiration date shown on the outside package label.
8. Handle the intraocular lens carefully. Rough handling or excessive handling may damage the lens.
9. The surgeon must be aware of the risk of opacification of the intraocular lens, which may necessitate lens removal. Although LENSTEC hydrophilic intraocular lenses have a satisfactory history regarding lens opacification, there is a history of lens opacification with lenses from other manufacturers. Most, if not all, of these types of cases required explanation. The material used by LENSTEC has not had any confirmed reports of 'Adverse Events' due to material discoloration, opacification and/or other material related deficiencies, which have caused postoperative patient problems.
10. All cases of lens removal must be reported to LENSTEC.
11. Medical facilities utilizing this IOL, and its accessories (if any) must ensure proper disposal as medical waste.

#### **HOW SUPPLIED**

The LENSTEC series of posterior chamber intraocular lenses are autoclave sterilized in a pre-loaded injector, contained within a double thermoform tray system. The contents of the outer tray are sterile unless the package is damaged or opened.

#### **INSTRUCTIONS FOR IMPLANTATION: SOFTEC PC IOL**

##### Calculation of Lens Power

It is recommended that the surgeon uses a power calculation method with which they are comfortable. In general, the power of the lens for each patient can be calculated from the keratometry measurements and axial length of the eye according to formulas in published literature. Additional reference to this topic can be found at [http://www.doctor-hill.com/iol-master/lens\\_constants.html](http://www.doctor-hill.com/iol-master/lens_constants.html)

NOTE: The 'A' Constant and ACD values printed on the outside of the package are estimates only. It is recommended that the surgeon determine his/her own values based on their individual clinical experience.

Pre-Surgical Preparation

- a. Determine the lens power from IOL Refractive Calculation Equation-Holladay or SRK/T.
- b. Determine the Expected Post-operative Target Refraction (SE)

**SURGICAL TECHNIQUE**

1. Peel the top of the small PLI Tip container (G) and carefully drop the PLI tip into the sterile field.
2. Open the outer tray of the PLI (F) and carefully drop its sterile inner tray into the sterile field.
3. Carefully peel the foil lid off the PLI inner tray and pour some saline out.
4. Remove the PLI from the tray and rinse with fresh saline – the red button (D) should be at the bottom.
5. Insert the PLI tip (C) into slot (B) at the front of the PLI.
6. Remove the red button from the PLI.
7. Press the clear button (A) firmly until it clicks into place.
8. Using the plunger (E), advance the lens into the PLI tip until the white silicone cushion is visible.
9. Perform safety check: Examine the PLI tip under a microscope while rotating it 360° to ensure that the trailing haptic is not caught between the cushion and the inner wall of the PLI tip.
10. If it is caught, retract the plunger fully, then advance it again. Repeat the safety check if haptic is still caught, until free.
11. The lens is now ready to be injected.

**EXPLANATION**

Explantation procedures may vary depending on patient condition and circumstances. The surgeon is therefore advised to use an explantation method which he/she determines will provide the most favourable patient outcomes.

**QUALITATIVE & QUANTITATIVE DATA FOR THE HEMA PLI MODELS**

HEMA PLI: Hydroxyethyl methacrylate, 26% water content. The devices have been tested and proven safe in accordance with ISO 10993-3, ISO 10993-5, ISO 10993-6, ISO 10993-7, ISO 10993-10, ISO 10993-11 and ISO 11979-5. Contact Lenstec for further details.

**DETAILED DEVICE DESCRIPTION**

Construction: Single Piece  
Material: 26% Water Content HEMA (Hydroxyethyl methacrylate)  
Light transmittance: Refer to diagram (H)  
Index of Refraction: 1.460

The specifications for the LENSTEC SOFTEC series of posterior chamber intraocular lenses are as follows:

Optic Size: 5.75 mm  
Optic Type: Equiconvex  
Length: 12.00 mm  
Angulation: 0 degrees  
Construction: 1 Piece  
Position Holes: 0 Holes  
Optic Material: HEMA (26% water content)

**EXPIRATION DATE**

The expiration date on the lens package is the sterility expiration date. Do not use the IOL after the expiration date.

**RETURNS POLICY**

Contact your Lenstec representative regarding the return goods policy. Return the lens with full identification and the reason for the return. Label the return package as a biohazard.

**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**

Refer to EUDAMED.

#### **PATIENT REGISTRATION AND REPORTING**

A Patient Identification Card is included in the package. This is to be completed and given to the patient, together with instructions to keep the card as a permanent record to be shown to any eye practitioner the patient consults in future. Self-adhesive lens identification labels are provided for use on the Patient Identification Card and other clinical records.

Adverse events/complaints that may reasonably be regarded as lens-related and that were not previously expected in nature, severity, or degree of incidence should be reported to the relevant EU Competent Authority of the Member State and Lenstec at Airport Commercial Centre, Pilgrim Road, Christ Church, Barbados: Tel: +1 246-420-6795 • Fax: +1 246-420-6797; Email: [feedback@lenstec.com](mailto:feedback@lenstec.com) or contact your Lenstec representative.

#### **BIBLIOGRAPHY (I)**

1. Holladay, J. Quality of Vision: essential optics for the cataract and refractive surgeon. Pgs 27-35
2. DJ, Arthur S, et al. Effect of heparin surface modification in reducing silicone oil adherence to various intraocular lenses. Journal of Cataract & Refractive Surgery. October 2001, Vol 27 No 10, 1662-1669
3. Hayashi H, Hayashi K, et al. Quantitative comparison of posterior capsule opacification after polymethylmethacrylate, silicone and soft acrylic intraocular lens implantation. Arch. Ophthalmol. 1998; 116: 1579-82
4. Apple DJ, Federman JL, et al. Irreversible silicone oil adhesion to silicone intraocular lenses. A clinicopathologic analysis. Ophthalmology, June 1996. Vol 103, 1555-1561
5. Pandey SK, Apple DJ, et al. Posterior Capsule Opacification: A Review of the Aetiopathogenesis, Experimental and Clinical Studies and Factors for Prevention. Indian J Ophthalmol 2004;52:99-112
6. Apple DJ, Kleinmann G, et al. A new classification of calcification of intraocular lenses. Ophthalmology. Jan 2008, Volume 115, Issue 1, Pages 73-79
7. ISO 11979-7 Ophthalmic Implants – Intraocular lenses – Part 7: Clinical investigations for intraocular lenses for the correction of aphakia; 2018.

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