http://www.lenstec.com
ISO 13485 Registered Company

LENSTEC (BARBADOS) INC., Airport Commercial Centre, Pilgrim Road, Christ Church BB17092, Barbados

+1 (246) 420-6795 +1 (246) 420-6797 lenstecbarbados@lenstec.com

CMI spol. s.r.o. Trenčianska 47,821 09 Bratislava, Slovakia

+421 (2) 52631441 +421 (2) 52922882 lenstec@cmi.sk


UK SALES: LENSTEC (BARBADOS) INC., Lenstec House, Unit 8, Mariner Court, Calder Park, Wakefield, WF4 3FL, UK

+44 (0) 1924 382 678 +44 (0) 1924 850 454 lenstecuk@lenstec.com

WORLD HEADQUARTERS: LENSTEC, INC., 1765 Commerce Ave N, Saint Petersburg, Florida 33716, U.S.A.

+1 (727) 571-2272 +1 (727) 571-1792 lenstec@lenstec.com

INSTRUCTIONS FOR USE • MODE D'EMPLOI • INSTRUCCIONES DE USO • ISTRUZIONI PER L'USO •
GEBRAUCHSANWEISUNG • INSTRUÇÕES DE UTILIZAÇÃO • NÁVOD K POUŽITÍ • POKYNY K POUŽITIU
• ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ • KULLANIM TALÍMATLARI • إرشادات الاستخدام • 使用説明 • GEBRUIKSAANWIJZING
• INSTRUCȚIUNI DE UTILIZARE

BASIC UDI

LR1300B: 08443690LR-1300B7K

EN	INSTRUCTIONS FOR USE	2
FR	MODE D'EMPLOI	4
ES	INSTRUCCIONES DE USO	6
IT	ISTRUZIONI PER L'USO	8
DE	GEBRAUCHSANWEISUNG	10
PT	INSTRUÇÕES DE UTILIZAÇÃO	12
CS	NÁVOD K POUŽITÍ	14
SK	POKYNY K POUŽITIU	16
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	18
TR	KULLANIM TALÍMATLARI	20
AR	إرشادات الاستخدام	22
ZH-Hans	使用説明	24
NL	GEBRUIKSAANWIJZING	26
RO	INSTRUCȚIUNI DE UTILIZARE	28

Sterilized using ethylene oxide

Do not reuse

Do not sterilize

Do not use if damaged

Caution, Consult accompanying documents

R_x only

CE
2797

PI08 Rev 11
12/09/2022

IMPORTANT NOTICE

This product is intended for use by a qualified ophthalmic surgeon in a sterile medical setting. It is highly recommended that the surgeon adhere to the recommendations, contraindications and warnings outlined in these instructions.

INTRODUCTION

Lenstec Capsular Tension Rings are non-optical devices and are intended to be used in combination with an intraocular lens for the surgical correction of the aphakia in adults 18 years and older. The rings are manufactured of medical grade Perspex CQ-UV Polymethylmethacrylate (PMMA).

INTENDED USE

Implantation of the Lenstec Capsular Tension Ring may be indicated in the following conditions:

- To stabilize the capsule in cases of defective or absent zonulae
- To prevent shrinkage/contraction of the capsular bag
- For circular expansion of the capsular bag
- To prevent/correct IOL decentration/tilt/dislocation
- To stabilize the capsule in cases of high myopia
- For stable operating conditions (e.g. after ocular trauma)

CONTRAINDICATIONS

These include the usual contraindications for cataract surgery and lens implantation, as well as the following operative complications in cataract surgery:

- Capsular bag defect
- Uncontrolled glaucoma
- Chronic severe uveitis
- Retina and optic nerve defects
- Diabetic retinopathy

CLINICAL BENEFITS / PERFORMANCE CHARACTERISTICS

The primary benefit of the capsular tension ring is to stabilize the capsular bag in cases of weak/absent zonules. The duration of the treatment effect is anticipated to be permanent.

RISKS

The potential risks of implanting the Lenstec Capsular Tension Rings are as follows: Endophthalmitis, Toxic anterior segment syndrome (TASS), Uveitis (including Iritis, Iridocyclitis, Vitritis), Anterior capsule fibrosis, IOL decentration/tilt, Elevated IOP, Concomitant surgery

Individuals in the following health categories may be at a higher risk of experiencing secondary IOL calcification (surface deposits on the lens):

- Diabetes
 - ◆ Associated:
 - * Retinal detachment
 - * Vitreous detachment
 - * Vitrectomy
 - * Diabetic retinopathy
 - * Diabetic maculopathy
- Glaucoma
- Underwent DMEK or DSEK (multiple procedures further increase the risk)
- Hypertension
- High cholesterol

WARNINGS

- Do not soak the capsular ring with any solution other than a sterile balanced salt solution or sterile normal saline.
- Do not use the capsular ring if the sterile pouch has been opened or damaged
- Do not re-use the ring
- Do not re-sterilize the capsular ring. Rings requiring re-sterilization should be returned to the manufacturer.
- Do not use the ring after the expiration date shown on the outside of the packaging. After this date, Lenstec cannot guarantee that the performance of the ring will remain unchanged.
- Inspect the ring prior to use
- Handle the ring with care. Rough or excessive handling may damage the ring.
- A high level of surgical skill is required for capsular tension ring implantation. A surgeon should have observed and/or assisted in numerous ring implant procedures and successfully completed one or more courses on cataract implants prior to attempting to implant the capsular ring.

- LENSTEC does not provide cleaning/sterilization instructions. An improperly cleaned and/or sterilized CTR can cause significant damage to a patient's vision, due in part to cross contamination induced infection.
- Once removed from its original packaging, the CTR can lose traceability. In the event a CTR is re-used, it is unlikely the user will know the correct expiry date, serial number or dioptic power.
- Medical facilities utilizing this device and its accessories (if any) must ensure proper disposal as medical waste.

DIRECTIONS FOR USE

Lenstec Inc. Capsular Tension Rings are ethylene oxide sterilized in a lens tray contained within two sealed Tyvek sterilizable peel pouches. The contents of the pouches are sterile unless the package has been opened or damaged.

Once the tray has been removed from the pouch, firmly hold the capsular tray in the palm of your hand and slide the lid open with your thumb. Inspect the capsular ring for debris and damage. Grasp the capsular ring along its circumference and rinse in a balanced saline solution prior to implantation.

The package box contains peelable labels which display the capsular ring serial number and model number. These labels are designed to be affixed to the patient's hospital chart and the surgeon's chart. One of these labels should be affixed to the patient's identification card contained in this package and given to the patient as a permanent record of their capsular ring implant.

IMPLANTATION

1. Make the incision for phacoemulsification. The size of this incision is purely the surgeon's preference based on the injector/cartridge system used.
2. Gently grasp the capsular tension ring at the side (closest to one of the ends), using a pair of non-toothed forceps, and remove it from the case.
3. Carefully introduce it into the capsular bag, and gently adjust the ring until the required position within the eye is obtained.
4. Upon alignment, implant the IOL (using the instructions for that device), carefully remove all viscoelastic, perform diligent polishing of the capsular bag, and ensure that the lens is positioned correctly before closing the wound.
5. Possible reasons for difficulty in implanting the capsular tension ring include:
 - a. Torn IOL
 - b. Defective/Incorrect injector/cartridge system
 - c. Capsular bag damage
 - d. Inadequate lubrication (e.g. viscoelastic)

EXPLANATION

Explanation procedures may vary depending on patient condition and circumstances. The surgeon is therefore advised to use an explanation method which he/she determines will provide the most favourable patient outcomes.

QUALITATIVE & QUANTITATIVE DATA FOR THE PMMA CTRs

PMMA (clear and blue): Medical Grade Perspex CQ-UV Polymethylmethacrylate. The devices have been tested and proven safe in accordance with ISO 10993-3 and ISO 11979-5. Contact Lenstec for further details.

EXPIRATION DATE

The expiration date on the lens package is the sterility expiration date.

RETURNS POLICY

Contact your Lenstec representative regarding the return goods policy.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The SSCP is available in the European database on medical devices (EUDAMED) where it is linked to the BASIC UDI. The EUDAMED website can be accessed at <https://ec.europa.eu/tools/eudamed>.

PATIENT REGISTRATION AND REPORTING

A Patient Identification Card is included in the package. This is to be completed and given to the patient, together with instructions to keep the card as a permanent record to be shown to any eye practitioner the patient consults in future. Self-adhesive lens identification labels are provided for use on the Patient Identification Card and other clinical records.

Adverse events/complaints that may reasonably be regarded as device-related and that were not previously expected in nature, severity, or degree of incidence should be reported to the relevant EU Competent Authority of the Member State and Lenstec at Airport Commercial Centre, Pilgrim Road, Christ Church, BB17092, Barbados: Tel: +1 246-420-6795 • Fax: +1 246-420-6797; Email: feedback@lenstec.com or contact your Lenstec representative.

IMPLANT CARD

- All fields present on the ID card must be completed by the healthcare institution/provider.
- It is the responsibility of the healthcare institution/provider to attach the label sticker with the etched 'LENSTEC' logo to the reverse (unprinted) face of the patient ID card, and provide this to the patient as a record of their implant.

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