

Ask About Younger Vision

*Near,
Far,
clear*





The **Tetraflex**TM
Presbyopic IOL
www.TheTetraflex.com



Younger Vision Available with


The **Tetraflex** Presbyopic IOL

- 1.** Your eye works much like a camera. Light enters through the cornea and passes through your natural lens, which bends these light rays so they focus on your retina, which transmits the image to your brain. Your natural lens subtly changes shape to clearly focus the image. The clarity and sharpness of the image on your retina determines how well you see.



The Tetraflex IOL is a new style of lens, designed to mimic the action of your natural lens and move subtly back and forth. This allows the image to be focused on your retina, giving you clear near, mid and distant vision once again, without glasses or contact lenses.

- 3.** With The Tetraflex IOL you will enjoy:
- reduced need for reading, driving or computer glasses
 - full peripheral (side to side) vision
 - a permanent solution to your vision problems
 - you can never develop cataracts with the removal of your natural lens



See Near, Far and Clear

www.TheTetraflex.com