

**LENSTEC (BARBADOS) INC.**, Airport Commercial Centre, Pilgrim Road, Christ Church BB17092, Barbados  
 +1 (246) 420-6795 +1 (246) 420-6797 lenstecbarbados@lenstec.com


**LENSTEC (SLOVAKIA) S.R.O.**, CMI spol. s.r.o. Trenčianska 47,821 09 Bratislava, Slovakia  
 +421 (2) 52631441 +421 (2) 52922882 lenstec@cmi.sk

**UK SALES: LENSTEC (BARBADOS) INC.**, Lenstec House, Unit 8, Mariner Court, Calder Park, Wakefield, WF4 3FL, UK  
 +44 (0) 1924 382 678 +44 (0) 1924 350 454 lensteouk@lenstec.com

**WORLD HEADQUARTERS: LENSTEC, INC.**, 1765 Commerce Ave N, Saint Petersburg, Florida 33716, U.S.A.  
 +1 (727) 571-2272 +1 (727) 571-1792 lenstec@lenstec.com

http://www.lenstec.com  
 ISO 13485 Registered Company

## Softec HEMA



CE  
2797



PK111 Rev 13  
TR/AR

21/07/2023

### BASIC UDI

Softec I: 08443690SOFTECIWL  
 Softec HD: 08443690SOFTECHD8G



Sterilized using steam



Do not reuse



Do not resterilize



Do not use if damaged






Caution. Consult accompanying documents

Rx only

<sup>1</sup> Holladay, J. Quality of Vision: essential optics for the cataract and refractive surgeon. Pgs 27-35

<sup>2</sup> Pandey SK, Apple DJ, et al. Posterior Capsule Opacification: A Review of the Aetiopathogenesis, Experimental and Clinical Studies and Factors for Prevention. Indian J Ophthalmol 2004;52:99-112

<sup>3</sup> Apple DJ, Kleinmann G, et al. A new classification of calcification of intraocular lenses. Ophthalmology, Jan 2008, Volume 115, Issue 1, Pages 73-79

MODEL	POWERS (D)	TOLERANCES (D)	INCREMENT (D)	ACCESSORIES
 <b>SOFTEC I™</b>	-5.0 to +10.0 +10.5 to +30.0 +31.0 to +36.0	$\pm 0.25$ $\pm 0.25$ $\pm 0.5$	1.0 0.5 1.0	<b>19011S Reusable Injector</b> Disposable Injector LC16i LC16, Cart45S -5.0D to +26.0D LC24 -5.0D to +36.0D 
 <b>SOFTEC HD™</b>	+5.0 to +10.0 +10.5 to +14.5 +15.0 to +25.0 +25.5 to +30.0 +31.0 to +36.0	$\pm 0.25$ $\pm 0.25$ $\pm 0.11$ $\pm 0.25$ $\pm 0.5$	1.0 0.5 0.25 0.5 1.0	<b>19011S Reusable Injector</b> Disposable Injector LC16i LC16, Cart45S +5.0D to +26.0D LC24 +5.0D to +36.0D 